

# Elbow CPM

## Artromot E2 Compact Elbow



### TECHNICAL SPECIFICATIONS

#### Range of Motion (ROM)

Extension/Flexion	-5 - 140°
Pronation/Supination (Manual Setup)	-90 - 90°

#### Functions

Motor Control	On/Off
Pauses	0 - 30 sec.
Speed	0 - 100%
Timer	1 - 300 min.
New Patient	
Warm-up Mode	50% of Programmed ROM
Synchronization	On/Off
Reverse-on-load	1 - 25
Therapy Run Time	h : m
Display Contrast	

#### Weight and Dimensions

Weight	35.2 lbs (16 kg)
Dimensions	30.7 x 40.9 x 23.6 x 31.5 in

#### SAFETY STANDARDS

CSA - 601-1-M90  
UL - 2601 - 1

**Lantz**  
Medical

(866) 236-8889