

Elbow CPM

Kinetec 6080 Elbow CPM



Anatomically correct passive motion of the elbow joint with or without synchronised pronation/supination of the forearm. It is convenient and easy to use in a bed or chair. Post-operative treatment using this device provides an effective deterrent to stiffness and loss of range of motion.

Indications

For the treatment of injuries and diseases of the elbow that limit movement of the joint.

Indications include the following:

- Open reduction and internal fixation of intra-articular fractures
- Prosthetic elbow joint replacement
- Arthroplasty
- Metaphyseal fractures in the elbow area with ORIF
- Arthrolysis for post-traumatic stiffness with limitation of elbow joint motion

Range of Motion

- Flexion: 135°
- Extension: 0°
- Pronation/Supination: 90°
- Speed (cycle of 0°-135°): 3 min 50 sec to 1 min

Specification

- Patient height: 140 to 190cm
- Weight: 21kg
- Dimensions: 72 x 65 x 130cm
- Voltage: 110/240 volts 50/60Hz

Lantz
Medical

(866) 236-8889